

The Case for Autofocusing Telescopes for Low Vision Bioptic Drivers

April 2025

The Benefits of Bioptic Telescopes for Low Vision Drivers

Bioptic telescopes are small optical devices mounted toward the top of regular eyeglasses. With a slight head tilt, the wearer can switch between normal vision and a magnified view, allowing them to see distant objects more clearly. These devices can help eligible individuals with visual impairments qualify for a restricted driver's license by improving their ability to read road signs, signals, and spot potential hazards. However, because bioptics prescribed for driving are intended only for distance viewing, they do not help drivers see dashboard displays.

Outdated Bioptic Driving Regulations

While most U.S. states will license individuals to use bioptic telescopes, regulations vary by state, and include differences in the design and magnification power of the bioptic, as well as driving restrictions such as speed and permitted times of day. Because drivers are not expected to adjust the device while operating a vehicle, many states require either a fixed-focus bioptic or prohibit adjusting a focusable one while driving. This limits the device's usefulness for other tasks and makes it difficult for users to view the dashboard. These regulations were established before autofocusing bioptics became available, and therefore do not reflect the new technology nor the drivers need to monitor dashboard displays.

The Challenge of Modern Car Displays

Modern vehicles feature complex dashboards with text, symbols, and dynamic digital displays. These screens provide critical information, such as speed, fuel levels, GPS navigation, and alerts. However, for bioptic drivers with low vision, these digital displays are difficult—if not impossible—to read.

The Autofocusing Advantage: Enhancing Safety and Efficiency

Autofocusing bioptics automatically adjust to the appropriate distance, allowing the driver to see both distant objects (like the road) and near objects (such as the car's dashboard) clearly without needing to manually adjust the device. Autofocusing bioptics offer a safer and more practical solution for bioptic driving than traditional fixed-focus devices, as they provide the full range of clear vision needed to support all the driver's vision needs.

Conclusion

While traditional bioptic telescopes are essential for enabling individuals with low vision to be eligible to receive a driver's license, they fall short in enabling the full range of clear vision that car dashboards and displays require. Autofocusing bioptics now offer the most compelling option to support the combined distance and near viewing needs of today's low vision driver. Their ability to automatically adjust their focus between the road and the dash ensures drivers can easily read both the road and their car's digital interfaces, enhancing both safety and efficiency. As vehicles continue to integrate advanced technology, autofocusing bioptics are becoming a compelling option for low vision drivers, allowing them to drive more confidently and safely with bioptic devices.