



VES[®] FALCON Autofocus Bioptic

5.5X Power
Now Available
10.0°
Field of View

The next generation of Ocutech's pioneering Autofocusing Bioptic Telescope, proven effective in NIH-funded clinical studies.

- The world's only autofocusing bioptic telescope
- Designed specifically for the visually impaired
- Provides the most natural, hands-free magnified vision
- Near, middle or far, images are immediately clear – just like normal vision
- 4x magnification and focusing range as close as 13 inches (0.33m)
- Crisp, bright, high-contrast wide-field Keplerian optics provide 12.5 degree field of view
- Operates 8 hours on a single charge
- Comfortable to wear all day long
- Ideal for seeing signs, faces, church, theater, TV, shopping, traveling, school, work, board games, card playing, computers and perhaps even driving*



CE

** Driving where legal and properly licensed*

TECHNICAL INFORMATION:

Optics: Monocular Keplerian 4x, 12.5 degree field of view
Interchangeable for right or left eye

Focus range: 13" (0.33m) to infinity

Eyepiece corrections: Adjustable eyepiece corrects for +/- 8D spherical refractive error
Eyepiece corrections indicated for Rx > +/-8D sphere or 3D cylinder

Frame mounting: Bridge mounting system
Completely adjustable for PD and viewing angle

Battery: Operates more than 8 hours on a full charge
Rechargeable USB Battery included
USB-C Battery Cable

Weight: 3.2oz (90 grams) including frame (excluding carrier lenses)

Frame options: Appropriate for all Ocutech metal frame styles

Carrier lenses: Prescribe the patient's habitual distance prescription including bifocals if normally worn

Accessories: Includes case, head-strap, cleaning cloth, fitting tools.
Slip-behind filters available; Clip-on filter available

Warranty: One year limited warranty for defects in manufacture and workmanship

Suggested patient price: \$3900 USD, excluding carrier lenses, customization charges or professional fees



Eyepiece is easily adjustable for refractive errors