

# OCUTECH® VES® *FALCON* Autofocus Biopic



**The next generation of Ocutech's pioneering Autofocusing Biopic Telescope, proven effective in NIH-funded clinical studies.**



- The world's only autofocusing biopic telescope
- Designed specifically for the visually impaired
- Provides the most natural, hands-free magnified vision
- Near, middle or far, images are immediately clear – just like normal vision
- Crisp, bright, hi-contrast Keplerian optics provide the widest fields of view possible
- 3x, 4x and 5.5x power options
- Operates 8 hours on a single charge
- Comfortable to wear all day long
- Choose Silver or Black chassis color
- Ideal for seeing signs, faces, church, theater, TV, shopping, traveling, school, work, board games, card playing, computers and perhaps even driving\*

*\* Driving where legal and properly licensed*

## TECHNICAL INFORMATION:

### **OPTICS:**

Monocular Keplerian telescope  
Interchangeable for right or left eye

### **POWERS, FIELD OF VIEW, FOCUS RANGE:**

3x	15° field of view	focuses to 13" (0.33m)
4x	12.5° field of view	focuses to 13" (0.33m)
5.5x	9.5° field of view	focuses to 15" (0.38m)

### **EYEPIECE CORRECTIONS:**

Adjustable eyepiece corrects for  
+/- 8D spherical refractive error

Eyepiece corrections indicated for  
Rx > +/-8D sphere or 3D cylinder

### **FRAME MOUNTING:**

Bridge mounting system. Completely adjustable  
for PD and viewing angle

### **BATTERY:**

Operates 8 hours on a full charge  
Rechargeable USB battery included

### **WEIGHT:**

3.2oz (90 grams) including frame  
(excluding carrier lenses)

### **FRAME OPTIONS:**

Appropriate for all Ocutech  
metal frame styles

### **CARRIER LENSES:**

Prescribe the patient's habitual  
distance prescription including  
bifocals if normally worn

### **ACCESSORIES:**

Includes case, head-strap,  
cleaning cloth, fitting tools.  
Slip-behind filters available;  
Clip-on filter available

### **WARRANTY:**

One year limited warranty for defects in manufacture  
and workmanship, 90 days for battery and cable.



*Eyepiece is easily adjustable for refractive errors*

**OCUTECH®**  
Vision Enhancing Systems

109 Conner Drive, Suite 2105 Chapel Hill, North Carolina 27514  
919 967-6460 • 800-326-6460 • 919 967-8146 Fax • www.ocutech.com