

OCUTECH®

VES® EXPLORER

*Small in size, big in optics,
and a pleasure to wear!*

FEATURES:

- *Wide field-of-view Keplerian optics!*
- *Available in 3x and 4x powers*
- *Hi-tech design and 5 color choices!*
- *Focusable for any distance*
- *Lightweight and comfortable (less than an ounce!)*
- *Great for kids and adults*
- *Fits all Ocutech metal frames including pediatric sizes*

The VES®-Explorer includes all of the great features of Ocutech's VES Keplerian telescope systems in this small and modern design! Appealing to both children and adults, its light weight makes the Explorer comfortable to be worn all day long!

Available in 3x and 4x powers, the VES Explorer offers the important combination of quality optics, a wide field of view and ease of use at any working distance.

Prescribed for individuals with 20/200 or better vision, the VES-Explorer comes complete with Ocutech's custom designed metal frames available in a range of styles, colors and sizes (even pediatric). The user's normal distance prescription (and even bifocals) can also be incorporated. Ocutech's innovative frame mounting system makes fitting, demonstrating and prescribing high quality bioptic telescope systems easy, convenient and risk-free! Even fine-tuning of telescope positioning after dispensing and future prescription lens changes are simple to make.

OCUTECH®
Vision Enhancing Systems

109 Conner Drive, Suite 2105 Chapel Hill, North Carolina 27514 • 919 967-6460 • 800-326-6460 • 919 967-8146 Fax • www.ocutech.com

© OCUTECH 2019



***Ocutech's online webinars will show you
just how easy it is to fit and prescribe
Ocutech's Bioptic Telescopes.***

OCUTECH®

VES® EXPLORER



TECH SPECS *Ocutech VES®-Explorer*

Power and field of view	3x, 14.0 deg	4x, 12.5 deg
Optical Design	Keplerian, monocular	
Telescope Chassis colors	Black, Silver, Blue, Red, Pink	
Interchangeable for right or left eye	Yes	
Focusing method	Central focusing ring	
Closest near focus (at neutral)	3x: 7 inches	4x: 9 inches
Weight (telescope alone)	0.7 oz (20 grams)	
Eyeglass frame options	Ocutech K, Unisex and Sleek frames only	
Mounting method	Adjustable bridge assembly	
Refractive error correction method	Internal focusing +12D/-12D	
Eyepiece corrections (suggestions)	Spheres above +/- 12 Cylinders indicated above 3D	
Adjustable after assembly	Yes	
Accessories	Slip-behind filters	
Warranty	2 years against defects	

OCUTECH®

Vision Enhancing Systems

109 Conner Drive, Suite 2105 Chapel Hill, North Carolina 27514 • 919 967-6460 • 800-326-6460 • 919 967-8146 Fax • www.ocutech.com

© OCUTECH 2019